

RESPIRONICS®



PEOPLE FIND COMFORT IN SIMPLE  
PLEASURES. THAT'S WHY WE MADE  
SLEEPEASY® AS SIMPLE AS POSSIBLE.



Envisioning tomorrow. Improving today.<sup>SM</sup>

**SleepEasy**®  
Sleep Therapy System

## SLEEPEASY SLEEP THERAPY SYSTEM: DESIGNED FOR TODAY'S SLEEP THERAPY BUSINESS.

Like all aspects of healthcare, the sleep therapy business is changing rapidly. But one thing hasn't changed: You need solutions that improve the health of your patients—and the health of your business. The SleepEasy sleep therapy system was designed to help you achieve both. It's a simple, reliable value-priced CPAP that you can trust—and that patients will use.



### VALUE: NOT JUST VALUE PRICED

Performance, portability and value are all within reach thanks to the combination of features, appealing design and price that you'll find only in SleepEasy.

### QUALITY YOU CAN COUNT ON

SleepEasy combines proven technologies and designs that have been Respironics standards for more than 30 years—we just made them more affordable.

### EASY TO USE

A permanently attached humidifier and clear, tactile control buttons make this sleep therapy system easy to set up and easier to live with.

### PRODUCT SPECIFICATIONS

#### SLEEPEASY

Pressure	4 to 20 cm H <sub>2</sub> O
Ramp Time	0 to 45 min (5-min increments)
Starting Ramp Pressure	4 to CPAP; patient-adjustable
Dimensions	25.4 x 19.1 x 14.0 cm
Weight	Approx. 1.6 kg
Filters	Pollen and optional ultra-fine
Device Set-Up	LED/Keypad
Data Storage Capacity (minimum)	Therapy meter
Data Storage Content	Total therapy hours
Total blower hours	Usage sessions > 4 hrs
Compliance Monitoring	Breathing detection
Altitude Compensation	Low, medium or high
Electrical Requirements	100–240 VAC, 50/60 Hz
Humidification	Fully integrated; patient-adjustable, 5 settings
Water Capacity	350 ml at recommended water level
Warranty	2 years
Misc.	Fine pressure calibration, lighted LEDs

### ORDERING INFORMATION

SLEEPEASY	PART NUMBER
SleepEasy (includes humidifier), International	1043238
SleepEasy (includes humidifier), Canada	1043239
ACCESSORIES	PART NUMBER
SleepEasy Replacement Chamber Kit	1048988
SleepEasy Carrying Case	1043517
Performance Tubing	1032907
Pollen Filter Kit	1029330
Ultra-Fine Filter Kit	1035442
Soft Valve Kit	1035323
Encore® Pro Software	1009583

**RESPIRONICS®**

Envisioning tomorrow. Improving today.™

Customer Service: +1-724-387-4000

Respironics Europe, Africa, Middle East: +33-1-47-52-30-00

Respironics Germany: +49-81-52-930-60

Respironics France: +33-2-51-89-60-00

Respironics Italy: +39-03-62-63-43-1

Respironics UK: +44-800-1300-845

Respironics Switzerland: +41-62-745-1750

Respironics Asia Pacific: +852-3194-2280

Respironics Australia: +61-2-9666-4444

CAUTION: U.S. federal law restricts these devices to sale by or on the order of a physician.

Respironics, SleepEasy and Encore are registered trademarks of Respironics, Inc. and its affiliates.

©2007 Respironics, Inc. and its affiliates. All rights reserved.