



## THE COMFORT YOUR PATIENTS NEED. THE VALUE YOU DESIRE.

Introducing ComfortFusion™. A value and a fit that are easy to get comfortable with. Our ComfortFusion mask has everything your patients want – a reliable seal, comfortable headgear, an easy fit and a replaceable cushion. Our FitPack™ comes with both a small and a medium cushion, enabling you to fit more than 90% of your patients – making this great mask a great value for you. It's the perfect blend of comfort and value – the ComfortFusion mask from Respiration.

ComfortFusion™



# ComfortFusion™

**Comfort and value come together in one dependable mask to benefit both you and your patients.**

The ComfortFusion mask, the newest member of the Respirationics Comfort Series™ of masks, combines our years of mask design experience with our understanding and appreciation of your business needs. Your patients will love the superior seal that results from the pairing of the soft, replaceable cushion and the StabilitySelector,® along with the premium soft forehead pad and the Fusion Headgear.™ You'll love that the ComfortFusion mask comes with two cushion sizes in our FitPack™— cushions that will fit more than 90%\* of your patients. Our DuoPack™ replacement kit comes with a mask (with cushion), plus an additional cushion of the same size.

The ComfortFusion mask has a familiar design that makes it easy to use, fit and achieve a robust seal. The result is a mask that's more comfortable for patients, which can lead to longer-term compliance. Overall, the ComfortFusion mask proves that comfort for your patients and value for your business belong together. And it's another example of how the Comfort Series family of masks from Respirationics helps fit 100% of your patients 100% of the time.

*\*Based on the Respirationics anthropometric database.*

## ORDERING INFORMATION

ITEM	PART NUMBER
ComfortFusion with Headgear FitPack Size S and M, Intl	104838
ComfortFusion with Headgear FitPack™ Size S and M, Intl	1042211
ComfortFusion with Headgear Size S, Intl	1042212
ComfortFusion with Headgear Size M	1040770
ComfortFusion without Headgear FitPack Size S and M	1040836
ComfortFusion with Headgear Size S - NO EXH, Intl	1044215
ComfortFusion with Headgear Size M - NO EXH, Intl	1044216
ComfortFusion with Headgear FitPack™ Size S and M, JPN	1040837
ComfortFusion with Headgear Size S, JPN	1042209
ComfortFusion with Headgear Size M, JPN	1042210
ComfortFusion with Headgear Size S - NO EXH, JPN	1040839
ComfortFusion with Headgear Size M - NO EXH, JPN	1040772
ComfortFusion with Headgear DuoPack Size S/S, Intl	1044213
ComfortFusion with Headgear DuoPack Size M/M, Intl	1042214

## REPLACEMENT PARTS


ITEM	PART NUMBER
ComfortFusion Cushion and Ring RP Kit Size S	1040840
ComfortFusion Cushion and Ring RP Kit Size M	1040841
Premium Silicone Forehead Pad 4 Per Pack RP Kit; Blue	1040774
Premium Silicone Forehead Pad with Support RP Kit 1 per pack; Blue	1040773
Fusion Headgear™ RP Kit only	1040844
ComfortFusion Sizing Gauge - Single	1041338
Premium Silicone FH pad RP Kit; Blue	1041960



The StabilitySelector raises and lowers the mask for optimized mask angle, seal and fit.

Two cushions (small and medium) fit more than 90% of patients and can be removed easily for cleaning and replacement.

The ergonomic ball-and-socket Quick Clips™ allow for free movement and easy mask removal.



The premium soft forehead pad provides superb patient comfort.

The soft, replaceable cushion provides a robust seal and stable fit.



The Fusion Headgear provides reliable comfort, fit and stability.

## RESPIRONICS®

Respirationics, Inc. : +1-724-387-4000  
 Respirationics Europe, Africa, Middle East : +33-1-47-52-30-00  
 Respirationics Germany : +49-8152-93-06-0  
 Respirationics France : +33-2-51-89-36-00  
 Respirationics Italy : +39-03-62-63-43-1  
 Respirationics UK : +44-800-1300-845  
 Respirationics Switzerland: +41-62-745-1750  
 Respirationics Asia Pacific: +852-3194-2280  
 Respirationics Australia: +61-2-9666-4444  
[www.respirationics.com](http://www.respirationics.com)

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician.

Respirationics and StabilitySelector are registered trademarks, and ComfortFusion, Comfort Series, DuoPack, FitPack, Fusion Headgear and Quick Clip are trademarks, of Respirationics, Inc. and its affiliates. Patents Pending.

©2007 Respirationics, Inc. and its affiliates. All rights reserved.

Broudy CY 10/28/07 MCI 4101361 PN 1002739