



Gloves wish
they fit this well.

The **ComfortGel™** Nasal Mask brings a new level of fit and comfort to the Comfort Series.™ An exciting new approach to fully customizable gel masks combines premium comfort features with a host of new innovations. These changes help the ComfortGel fit right the first time. And that helps the Comfort Series fit 100% of your patients 100% of the time.

t h e c o m f o r t s e r i e s ™

ComfortGel™

Comfortable and customizable, the latest addition to the Comfort Series™ of Masks fits right the first time.

The ComfortGel™ Nasal Mask makes it easier to fit patients right the first time – even hard-to-fit patients. At its core is an exceptional cushion that combines the customizable blue gel originally seen in the Profile™ Lite with a proven dual-layer cushion design found in other Comfort Series masks. This combination gives you the comfort, fit, stability and seal characteristics you've come to expect from Comfort Series masks.

This superior cushion is paired with the new StabilitySelector™ that tilts the mask to allow for an optimal seal. The forehead pad automatically adjusts to ensure proper forehead placement and angle. Combining the StabilitySelector and the dual material cushion with a silent exhalation port and new "ball-and-socket" headgear attachments adds up to a superior mask that can help your patients stick to their therapy.

Delivering more features in a smaller, more compact package makes the ComfortGel an easy choice. Our streamlined sizing system makes it easier to stock. And adding the ComfortGel to the Comfort Series makes it possible to fit 100% of your patients 100% of the time.

To learn more about the ComfortGel and all of the masks in the Comfort Series, visit the Comfort Series website:

www.comfortseries.respironics.com


ordering info

Item	Petite	Small	Med	Large
COMFORTGEL MASK WITH HEADGEAR	1009044	1009045	1009046	1009047
COMFORTGEL MASK WITH HEADGEAR, NO EXHALATION PORT	1010642	1010643	1010644	1010645

Custom-Sized Accessories	Petite	Small	Med	Large
GEL CUSHION	1009048	1009049	1009050	1009051
SILICONE COMFORT FLAP	1009052	1009053	1009054	1009109

One-Size-Fits-All

FOREHEAD SUPPORT ARM	1009721
SILICONE FOREHEAD PADS (4-PACK)	1009722
HEADGEAR SWIVEL CLIPS (2-PACK)	1009723
EXHALATION ELBOW	1010646
MASK SWIVEL	7041
DELUXE HEADGEAR (BLK-4-STRAP)	1002801



ComfortGel uses Respironics' one-size-fits-most **Deluxe Headgear**. This unique headgear with Velcro® tabs easily adjusts in eight places for a secure, custom fit. The Deluxe Headgear evenly distributes tension, adds stability and helps to prevent mask leaks. ComfortGel and Headgear disassemble easily for cleaning and replacement. The Headgear can be used with most Respironics' masks.



The ComfortGel's blue gel inner layer can be easily customized to give patients a truly personal fit.



Respironics International: +33-(0)1-55-60-19-80
 Respironics France: +33-(0)2-51-89-60-00
 Respironics Deutschland: +49-81-52-930-60
 Respironics Asia Pacific: +852-234-342-18
www.respironics.com



CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician. ComfortGel, Comfort Series, Comfort Flap, StabilitySelector and Profile Lite are trademarks and "People. Products. Programs." is a service mark of Respironics, Inc. and its affiliates. Velcro® is a trademark of Velcro Industries. Latex Free • Patent Pending

©2003 Respironics, Inc.

1012078 KW 03/20/03

